



# BIOCLEAN LEGACY™ Non-sterile Latex Gloves BLA2

## About

BioClean Legacy™ Latex Gloves provide maximum wearer comfort with exceptional flexibility. They feature a beaded cuff for strength and stability on the arm and a textured surface for enhanced grip.

## Specifications

**COMPATIBILITY:**  
ISO Class 4

**LENGTH:**  
300mm (12")

**MATERIAL:**  
Latex

**PROTEIN LEVEL:**  
Less than 50µg/gram

**SURFACE:**  
Textured


**SHAPE:**  
Ambidextrous

**COLOUR:**  
Natural

## Features

- Exceptional flexibility and comfort
- Processed to ensure ISO Class 4 compatibility
- Textured for improved grip
- Beaded cuff for strength and stability on the arm

POWDER-FREE  
TEXTURED  
BEADED CUFF  
UNBEATABLE COMFORT



Ms. G. Love™ represents BioClean's range of cleanroom gloves

## Quality Standards

- Manufactured in a facility holding ISO 9001: 2008 and 89/686/EEC Article 11b certifications
- Conforms to Category 3 - Complex Design Personal Protective Equipment PPE Directive 89/686/EEC
- Processed in an NEBB certified ISO Class 4 environment
- Complies with the requirements of European Glove Standards EN420:2003 and EN374-1, 2 & 3:2003
- Physical properties comply with European Medical Glove Standard EN455-2:2009

## Shelf Life & Storage

Five (5) years from date of manufacture.  
Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

## Physical Properties

CHARACTERISTICS	VALUE	TEST METHOD	PALM WIDTHS (mm) ± 5mm						
Freedom From Holes	AQL 1.5	EN 374-2:2003	SIZE	XSmall	Small	Medium	Large	XLarge	XXLarge
Minimum Length (mm/inches)	290/11.4	EN 420:2003	mm	77	83	95	106	112	118
Single Wall Thickness Minimum (mm)	Cuff: 0.11 Palm: 0.17 Finger Tip: 0.20	EN 455-2:2009	TEST METHOD				EN455-2:2009		
Force At Break During Shelf Life	≥ 9 N	EN 455-2:2009							

## Glove Cleanliness Characteristics

PARTICLES				TYPICAL PARTICLE COUNT				TEST METHOD
≥ 0.5µm(counts/cm <sup>2</sup> )				<1500				IENT-RP-CC005.3
TYPICAL EXTRACTABLE IONS								TEST METHOD
ANIONS	Fluoride	Chloride	Nitrite	Bromide	Nitrate	Phosphate	Sulphate	IENT-RP-CC005.3
(µg/cm <sup>2</sup> )	ND	0.97	NT	NT	0.67	0.07	ND	
CATIONS	Lithium	Sodium	Ammonium	Potassium	Calcium	Magnesium	Zinc	
(µg/cm <sup>2</sup> )	NT	0.02	NT	0.03	0.37	0.003	ND	
ND = Not Detected, NT = Not Tested								

## TO ORDER

### Re-order code

BLA200-XS  
BLA201-S  
BLA202-M  
BLA203-L  
BLA204-XL  
BLA205-XXL

### Size

XSmall  
Small  
Medium  
Large  
XLarge  
XXLarge

### Packing

100 pieces per sealed inner PE bag;  
one inner PE bag per sealed outer PE  
bag; 10 outer bags per lined carton  
(1000 pieces)

## CONTACT

### Europe (Headquarters)

Nitritex Ltd.  
Minton Enterprise Park,  
Oaks Drive,  
Newmarket, Suffolk,  
CB8 7YY,  
United Kingdom  
T: +44 (0) 1638 663338  
F: +44 (0) 1638 668890  
E: [europa@bioclean.com](mailto:europa@bioclean.com)

### Other Regions

Asia: [malaysia@bioclean.com](mailto:malaysia@bioclean.com)  
Americas: [americas@bioclean.com](mailto:americas@bioclean.com)



BioClean™ is a brand of Nitritex Ltd.

For more information, please visit [www.bioclean.com](http://www.bioclean.com)

BLA2/1603/PDS19